


REUSABLE FOREHEAD THERMOMETER



To order, contact:
sales@denverathletic.com



#9424 REUSABLE FOREHEAD THERMOMETER

- 1 Thermometer Per Case
- Provides Core Temperature Quickly And Easily
- Liquid Crystal Display Stops At Core Temperature For Quick Temperature Reading
- Use Indoor And At Room Temperature
- Avoid Taking Temperature In Direct Sunlight
- Reusable • Available May 1st
- Meets FDA Requirements

As Low As **\$1.75**

Reusable Forehead Thermometer Instructions:
USA Made FDA registered medical device. Use indoor and at room temperature (approximately 68-86 degrees F). Hold liquid crystal thermometer at both ends, pressing firmly on top forehead for 15 seconds. Watch the crystal display until it stops changing colors. Read thermometer while on forehead. Green is correct temperature. If you do not see green then blue or tan colors will appear. Always consult your physician if a fever exists. The reusable thermometer may be cleaned with a soft cloth and should be stored inside instruction packaging.
To be administered by adults only, keep out of reach of children.

